

# Atomic Layer Deposition

4D Labs system and recipes

David Dvorak

June 2016

## Cambridge NanoTech Fiji F200

From 4D Labs website:

The ALD system is equipped with a heated process chamber, remote plasma source, turbomolecular pump, and automated load lock for transferring substrates. The system is primarily intended for the deposition of a variety of thin films, including metals, oxides, and nitrides. Typical film thicknesses are monolayers up to tens of nanometers. Thicker films may be deposited but, since typical deposition rates are  $\sim 3\text{-}20\text{ nm/h}$ , deposition times will be several hours, at least.

### **Films available:**

Currently, the following films may be deposited in the ALD system; however, it has the capability to deposit many other materials, which may be added on request. For additional information, please contact the tool owner.

Oxides:  $\text{Al}_2\text{O}_3$ ,  $\text{HfO}_2$ ,  $\text{SiO}_2$ , &  $\text{ZnO}$ .

Nitrides:  $\text{AlN}$

Metals:  $\text{Pt}$

### **Key features:**

Highly conformable, well controlled layer by layer film deposition.

Multi-layer film stacks possible.

Aspect ratios of up to 450:1 in thermal mode or 20:1 in plasma mode are possible.

Substrates up to 200 mm in diameter and 6 mm thick can be accommodated.

Substrates may be heated up to  $400^\circ\text{C}$ .

Process chamber turbomolecular pump provides a base pressure  $\sim 2 \times 10^{-5}$  Torr.

Load lock chamber with automated transfer of substrates into the process chamber.


All depositions are performed from recipes.

Available plasma gases are  $\text{Ar}$ ,  $\text{O}_2$ ,  $\text{N}_2$ , &  $\text{H}_2$ . Any combination of gases may be used simultaneously, with the exception that  $\text{O}_2$  &  $\text{H}_2$  cannot flow simultaneously.






# $\text{Ga}_2\text{O}_3$ – New TMGa source Sept 2015





# $\text{Ga}_2\text{O}_3$ – Purge Headspace

- Old purge headspace recipe, data from May 2015





# $\text{Ga}_2\text{O}_3$ – Purge Headspace

- New purge headspace pressure data from Sept. 2015 (right)
- Purge headspace recipe on file, not sure if it's old or new (left)

Recipe: ALD5\_Purge\_Headspace

Status: Changed but not saved

	Instruction	#	Value	Units
0	APC		100	%
1	flow	0	40	sccm
2	flow	1	40	sccm
3	wait		10	s
4	pulse	5	1	s
5	wait		5	s
6	goto	4	20	cycles
7	flow	0	10	sccm
8	flow	1	30	sccm

MFC 0 Ar Carrier (sccm)  
10 9.93

MFC 1 Ar Plasma (sccm)  
30 29.4

MFC 2 N2 Plasma (sccm)  
0 1.81


MFC 3 O2 Plasma (sccm)  
0 1.78

MFC 4 H2 Plasma (sccm)  
0 1.39

MFC 5 No gas  
0 1.47

MFC 6 No gas  
0 1.51

MFC 7 No MFC  
0 0







# Ga<sub>2</sub>O<sub>3</sub> – Growth

- Same recipe used every time, only # of cycles is changed

Recipe: Ga<sub>2</sub>O<sub>3</sub>\_350C\_TMGa\_O2\_Plasma

Status: Changed but not saved


	Instruction	#	Value	Units
15	MFCvalve	3	1	
16	wait		20	s
17	pulse	5	0.015	s
18	wait		5	s
19	flow	3	25	sccm
20	wait		1	s
21	plasma		300	W
22	wait		10	s
23	plasma		0	W
24	flow	3	0	sccm
25	wait		5	s
26	goto	17	189	cycles
27	flow	0	10	sccm
28	flow	1	30	sccm






# $\text{SiO}_2$ – Growth & Purge Headspace

Growth




Purge Headspace






# $\text{Al}_2\text{O}_3$ – Growth & Purge Headspace

Growth




Purge Headspace



# Summary

	Old Ga <sub>2</sub> O <sub>3</sub>	Ga <sub>2</sub> O <sub>3</sub> (Sept 2015 ->)	Al <sub>2</sub> O <sub>3</sub>	SiO <sub>2</sub>
Purge Headspace – spike height	1.0	0.46-0.52	0.14	0.03
Purge Headspace – spike peak pressure	0.01	0.01	0.15	0.04
Growth– spike height	0.06	0.07	0.02	0.015
Growth – spike peak pressure	0.16	0.16	0.15	0.09





# ALD Runs:

Precursor	Date	Nominal Thickness (nm)	# of cycles	ALD ticket (hrs)	Cleanroom ticket (hrs)
TDMASi	April 8, 2015	8		5.58	1.67
TMGa	May 1, 2015	10	189	4.47	2.25
TDMASi	May 8, 2015	2.5	139	2.85	2
TMAI	June 26, 2015	8	90	2	1
TMGa	Sept. 8, 2015	8	94	2.33	1
TMGa	Sept. 16, 2015	8	93	2.32	1
TMGa	Dec. 21, 2015	40	465	5	3
TMAI	Jan. 13, 2016	40.5	450	3.17	1.25
TMAI	Jan. 18, 2016	8	90	1.5	0.5
TMAI	Jan. 21, 2016	2	22	1	0.5
TMAI	Jan. 21, 2016	4	45	1.33	0.5
TMGa	Feb. 2, 2016	20	233	2.25	0.75
TMAI	Feb. 23, 2016	2	22	1.5	0.5
TMAI	June 8, 2016	2	22	1.33	1

### Deposition Rates For 4D Recipes

Film	Recipe	Precursor #1		Precursor #2	T <sub>chuck</sub> °C	No. Cycles	Cycle Time (s)	Deposition	
		Material	T (°C)					nm/cycle	nm/h
AlN	AlN_200C_TMA_N2-H2-4D-20111123 *	TMA	RT	N <sub>2</sub> -H <sub>2</sub>	200	800	56.06	0.075	4.8
Al <sub>2</sub> O <sub>3</sub>	Al2O3_250C_TMA_H2O-4D-20111123	TMA	RT	H <sub>2</sub> O	250	1,000	20.12	0.089	15.9
HfO <sub>2</sub>	HfO2_200C_TDMAH_H2O_PRK_20111212	TDMAHf	75	H <sub>2</sub> O	200	1,000	20.3	0.108	19.1
Pt	Pt_270C_MeCpPtMe3_O2-4D-20111202	(MeCp)PtMe <sub>3</sub>	80	O <sub>2</sub>	270	350	26.0	0.043	5.9
Pt	Pt_200C_MeCpPtMe3_O2-4D-20111202	(MeCp)PtMe <sub>3</sub>	80	O <sub>2</sub>	200	300	36.0	0.036	3.6
Pt	Pt_150C_MeCpPtMe3_O2-4D-20111202	(MeCp)PtMe <sub>3</sub>	80	O <sub>2</sub>	150	300	66.0	0.036	2.0
Pt	Pt_90C_MeCpPtMe3_O2-4D-20111202	(MeCp)PtMe <sub>3</sub>	80	O <sub>2</sub>	90	300	96.0	0.036	1.4
SiO <sub>2</sub>	SiO2_250C_TDMASi_O2-4D-20111123	TDMASi	RT	O <sub>2</sub>	250	1,200	19.1	0.018	3.4
SiO <sub>2</sub>	Passivation-SiO2_250C_TDMASi_O2-4D-20111123	TDMASi	RT	O <sub>2</sub>	250	1,400	19.1	0.018	3.4
ZnO	ZnO_250C_DEZ_H2O-4D-20111123	DEZ	RT	H <sub>2</sub> O	250	800	20.12	0.130	23.3

\* Do not use this recipe. Films degrade over time.  
Deposition rate in regular face type have been measured at 4D Labs but rates in italics are based on Cambridge NanoTech standard recipes.